



Canadian Food  
Inspection Agency


Agence canadienne  
d'inspection des aliments

# Nanopore Basecalling and Read QC



**GRDI-AMR2  
Metagenomics Working Group**

**Marc-Olivier Duceppe  
April 22, 2024**

# Presentation Plan

1. Why Nanopore?
  2. How to basecall the reads in “super accuracy” mode?
  3. How to assess run / read quality?
- 

# The Nanopore Technology



- Basecalling is AI based (Deep Neural Network)
- Performed from 5-mers

<https://nanoporetech.com/platform/technology/basecalling>

# Two Families of Sequencers




- MinION Mk1B
- Flongle
- GridION
- **Cheaper**



- PromethION
- P2 Solo
- P2i
- **Higher yield and Q-score**

<https://store.nanoporetech.com/>

# Do You Believe in Nanopore?



## Pros


- \$1k / 87g
- **20bp - 100+kb**
- Native DNA / RNA
- Real-time sequencing
- Reusable flowcells
- Adaptive sequencing
- Growing community

## Cons

- Keeps changing (QA)
- **Error-prone**
- Hard to analyze
  - Computationally intensive (GPU)
  - Somewhat specialized tools



# Accuracy



- Can't quite resolve with accuracy long homopolymers stretches
- Leads to lots of INDELs
- Often long with more errors is better than short with less errors
- AMR SNP-based: ✓
- AMR functional gene: ✗


# Evolution of Basecallers

- Albacore: 2017 – 2018
- Guppy: 2018 – 2023
- Dorado: 2023 – present



Wick RR, Judd LM, Holt KE. Performance of neural network basecalling tools for Oxford Nanopore sequencing. *Genome Biol.* 2019 Jun 24;20(1):129. doi: 10.1186/s13059-019-1727-y. PMID: 31234903; PMCID: PMC6591954.

# Read Accuracy Keeps Improving



- In 2024
  - R10.4.1 v14
  - Dorado v0.6
  - Super accuracy mode
- Simplex: 99.5% (Q23)
- Duplex: >99.9% (Q30)
- **More Q15 – Q20**

<https://nanoporetech.com/platform/accuracy>




# To Guppy or not to Guppy



- Why you'd still want to use Guppy
  - Don't know how to run Dorado
  - You have an easier way to run Guppy
  - Part of your current SOP / certified workflow
  - Downstream analysis requirements
    - "Old" tools with "old" data (e.g. base modification)
  - Lower end GPU (< 8GB VRAM)
- Why you don't want to use Guppy anymore
  - Legacy (no more updates)
  - Slower
  - Less accurate
  - New chemistries won't be supported

# Basecalling Recommendations

- While sequencing: “Fast” mode
  - Enables more real-time run metrics
  - Do not overload host computer
- Post sequencing: “Super accuracy” mode
  - Higher Q-scores
  - Less “fail” reads
  - Better barcode assignment
  - Need computer with high-end GPU
  - Lose real-time benefit



# Guppy - Legacy

- Input: Fast5 and kits (“final\_summary”)
- Output: Fastq and “sequencing\_summary”
- [https://cdn.oxfordnanoportal.com/software/analysis/ont-guppy\\_6.5.7\\_linux64.tar.gz](https://cdn.oxfordnanoportal.com/software/analysis/ont-guppy_6.5.7_linux64.tar.gz)
- [https://github.com/duceppemo/basecall\\_nanopore](https://github.com/duceppemo/basecall_nanopore)



# Dorado – Two Flavours

## Standalone

- Input: Pod5 | Fast5
- Output: a single uBAM | SAM | Fastq
- Needs to be demultiplexed separately
- Fastq headers only contain “ReadID”
- <https://github.com/nanoporetech/dorado> (v0.6.0)



# Dorado – Two Flavours



## Basecall Server

- Input: Pod5 | Fast5
- Output: Fastq and “sequencing\_summary.txt”
- Demultiplexing built in
- Fastq headers with “all” info (like Guppy)
- [https://cdn.oxfordnanoportal.com/software/analysis/ont-dorado-server\\_7.3.9\\_linux64.tar.gz](https://cdn.oxfordnanoportal.com/software/analysis/ont-dorado-server_7.3.9_linux64.tar.gz)



# Nanopore read QC



- pycoQC v2 (<https://github.com/a-slide/pycoQC>)
  - Designed for Guppy
  - Can work with Dorado
    - Standalone: need to run additional commands to generate a compatible "seq\_summaray.txt" file
    - Server: as is





# Added Features

- pycoQC v3 (<https://github.com/duceppemo/pycoQC>)
  - Work in progress
  - Designed to accept all dorado output file types
  - Additional scripts:
    - Pod5\_to\_seq\_summary (no basecall information)
    - Bam\_to\_seq\_summary (slow)
    - **Fastq\_to\_seq\_summary (recommended)**
  - Starting from Fastq triggers additional GC plots
  - Display “All”, “Pass” or “Fail” reads





# Suboptimal Loading

- High quality library
- Not enough loaded
  - Ran out after 10h
- High quality library
- Run stopped too early



# Suboptimal Sample Normalization

- Uneven barcode distribution
  - Poor normalization
  - Adapter ligation problems
    - Wrong molarity
- Equal barcode distribution
  - MinKNOW now uses adaptive sequencing to help balance barcodes



# Run Diagnosis

- Very useful to find what went wrong
  - DNA extraction (contaminants, short DNA fragments, etc.)
  - Library prep (adapter-to-DNA ratio, sample normalization, etc.)
  - Sequencing (flowcell quality, amount loaded, overheating, script error, etc.)



# You Can Do It!

- Github repo with basic code to install and run Dorado basecall server
- [https://github.com/OLF-Bioinformatics/2024-04-22\\_GRDI\\_basecalling\\_presentation](https://github.com/OLF-Bioinformatics/2024-04-22_GRDI_basecalling_presentation)

# Nanopore Metagenomics Resources

- <https://github.com/duceppemo/mashID>
  - Quick species ID from metagenomics samples
  - Raw reads or assemblies
  - Pre-compiled DB available for bacteria
- <https://github.com/duceppemo/seqcounter>
  - Ultra-rapid preliminary AMR detection from raw reads

# Take-Home Message

- Switch to Dorado!
- Why basecall your runs using Dorado SUP?
  - Access higher accuracy
  - It keeps getting better (chemistry and software)
  - Better downstream analysis results
- DNA quality (and its size) matters
  - Garbage in = garbage out
  - The longer the better
  - Amplified DNA sequences better (number of reads and q-score)
- Bioinformaticians are not magicians
  - Bad runs need to be scraped (like for diagnostics)
- The truth is in the reads
  - Read QC will help you find what went wrong and areas of improvement



**Do it.**



# Thank You!

## Q & A